



# CATALOGUE



# Diagnostic A

## EXAM SET UP 7 INS (Recommended)

1. 1007 3EZ EXPLORERS TU17/23
2. 1018 3EZ EXPLORERS ODU11/12
3. 1104 1EZ CP12
4. 6517 FORCEP LONDON-COLLEGE
5. 6911 6EZ MIRROR HANDLE #6 C/S
6. 6951 #5 MIRROR
7. 7501 6 HOLD CASSETTE



### HANDLE OPTION


- \* 1EZ S/E
- \* 3EZ D/E
- \* 6EZ D/E OR S/E

Explorers  
Probes(Color)  
Explorer/Probes  
Mirror Handles

To detect and explore caries, abrasions, calculus, pocket characteristics, furcation and restoration. Instruments do not bend easily due to our outstanding tip strength.

It also gives dentists a delicate feeling to make an exact diagnosis of the patient's condition.


More moderate tips are usually used for disease detection while hard tips are used for removing calculus and excessive cement.



1007 TU17: Double explorer comprised #17 for gingival calculus detection and #23 for carv detection




To check the depth of periodontal pocket vertically or horizontally. Length of 6mm is the standard control point for diagnosing patient's periodontal symptom: (a) If the length of periodontal pocket is below 6mm, patient is at the mild condition. (b) If the length of periodontal pocket is more than 6mm, patient is at the critical condition.





To check the depth of periodontal pocket vertically or horizontally. Length of 6mm is the standard control point for diagnosing patient's periodontal symptom: (a) If the length of periodontal pocket is below 6mm, patient is at the mild condition. (b) If the length of periodontal pocket is more than 6mm, patient is at the critical condition.



Variety of mirror sizes, styles, and handles.

MIRROR HANDLE C/S TYPE



AVAILABLE  
 6950 #4 Mirror Head Front Surface  
 6951 #5 Mirror Head Front Surface  
 6951D #5 Double Mirror Head Front Surface



Variety of mirror sizes, styles, and handles.



6931 PARA ANTERIOR 5mm



6932 PARA POSTERIOR 5mm




6933 ELEVATOR C/S MIRROR



# Periodontal B

## PROPHY SET UP 15 INS (Recommended)

1. 1018 3EZ EXPLORERS OUD11/12
2. 1132 3EZ COLOR PROBE EXP23/CP12
3. 1201 6EZ SCALERS/CURETTES H6/7
4. 1208 6EZ SCALERS/CURETTES 204S
5. 1416 6EZ SCALERS/CURETTES COL13/14
6. 1601 6EZ STANDARD GRACEY GR1/2
7. 1604 6EZ STANDARD GRACEY GR7/8
8. 1606 6EZ STANDARD GRACEY GR11/12
9. 1607 6EZ STANDARD GRACEY GR13/14
10. 6517 FORCEP LONDON-COLLEGE
11. 6911 6EZ MIRROR HANDLE #6 C/S
12. 6950 #4 MIRROR C/S FRONT SURFACE
13. 7423 SHARPENING STONE ARKANSAS
14. 7424 SHARPENING STONE ARKANSAS EDGE
15. 7502N CASSETTE 16 HOLD DOUBLE DECK




Sickle/Jacquette Scalers  
Black Titaniums  
Sickle Scalers/Curettes  
Scaler Curettes  
Universal Curettes  
Gracey Curettes  
Debridement/Furcations  
Sharpening Stones

HANDLE OPTION

\* 6EZ D/E OR S/E



To remove calculus on the surface of gingival area.



To remove calculus on the surface of gingival area.



1201  
H6/7





BLACK TITANIUM INSTRUMENTS

Our black titanium instruments are enhanced in strength and durability up to 3 times more than the ordinary tips. It is highly hygienic as it protects the instruments from virus infection. It is easy to clean and requires less sharpening than other instruments.

The titanium provides the surface of the instrument to be a special non-stick design.




1601T  
GR1/2  
BLACK TITANIUM

BLACK TITANIUMS



1201T  
H6/7  
BLACK TITANIUM



1203T  
H5/33  
BLACK TITANIUM



1208T  
204S  
BLACK TITANIUM



1226T  
NV-2  
BLACK TITANIUM



1322T  
MC13/14S  
BLACK TITANIUM



1415T  
COL4R/4L  
BLACK TITANIUM



1416T  
COL13/14  
BLACK TITANIUM



1422T  
BH5/6  
BLACK TITANIUM



1601T  
GR1/2  
BLACK TITANIUM



1606T  
GR11/12  
BLACK TITANIUM



1607T  
GR13/14  
BLACK TITANIUM



1608T  
GR15/16  
BLACK TITANIUM



1609T  
GR17/18  
BLACK TITANIUM



1610T  
AF1/2  
BLACK TITANIUM



1614T  
AF11/12  
BLACK TITANIUM



1615T  
AF13/14  
BLACK TITANIUM



1620T  
SUB1/2  
BLACK TITANIUM



1621T  
SUB11/12  
BLACK TITANIUM



1622T  
SUB3/4  
BLACK TITANIUM


Designed for the removal of supragingival calculus.



1326  
ALSO AVAILABLE MC17/18 WITH POINTED BLADE


SCALER & CURETTES

To remove heavy calculus above or below the surface of gingival area.  
 We angulate and sharpen the blade of our instruments more, to greatly reduce treatment fatigue.





Curettes are types of scalers with which dentists can deliver more precise treatments. The end of the tips are designed with two cutting edges and rounded toes.



LANGER (UNIVERSAL)



PROPHY GRACEY  
SHANKS ARE A LITTLE SHORTER THAN A NORMAL GRACEY



1505 YG7/8 The end of tips is designed with two cutting edges and rounded toe. It is used on all surfaces


■BUNT: Bunting ■YG: Younger-Good ■BT: Bates ■R: Rule ■MAL: Mallery ■HU: Hutchinson

GRACEY CURETTES


For the removal of heavy subgingival calculus including deep scaling, root planning and periodontal debridement. It is highly recommended for dentists to be equipped with a full set of Standard Gracey.

STANDARD GRACEY

RIGID & SHORTER  
BLADE AVAILABLE



SUB-0 GRACEY



FACIAL / LINGUAL

AF-GRACEY (AVAILABLE SHORTER BLADE)



SUB-0 GRACEY  
SMALL BLADE








MESIAL



DISTAL

- \* O' Hehir curettes debridement : Removing residual calculus after ultrasonic scaling for deep pockets.
- \* Furcation curettes : Used for curetting the infected areas of periodontal disease.







- \* Surgery Knives : For removal of soft tissue.
- \* Surgery Files : For removal of granulatous.



\* Surgery Knives : For removal of soft tissue.

\* Surgery Files : For removal of granulomatous.



1811  
SOLT1/2



1812  
KOR1/2



1812M  
KOR1/2M



1812L  
KOR1/2L



1813  
KCK-3



1816D  
KM3/4

- \* Surgery Knives : For removal of soft tissue.
- \* Surgery Files : For removal of granulosomatous.



1818  
F9/10S



1819  
F11/12S



1820  
F1/2S



1821  
F3/4S



For removal of remaining tissue of infected gingival area after extraction.



KIRKLAND



SUGARMAN



KRAMER NEVINS CURETTES



PERIOTOME (SINGLE END AVAILABLE)



For removal of teeth with minimal tissue damage.  
Serrated blade easily cuts periodontal ligament.



1923T  
PT-2 TITANIUM TG  
A/V Regular



1924T  
PT-1 TITANIUM TG  
A/V Regular



1925T  
PT-3 TITANIUM Angel S/E



1926T  
PT-4 TITANIUM STR S/E



1927T  
PT-5 TITANIUM TG CVD S/E



1928  
MS-1 PERIOTOME



2043T  
PPAEL TITANIUM  
REGULAR AVAILABLE



2044T  
PPAELA TITANIUM  
REGULAR AVAILABLE

For scraping out and reshaping bone.  
For removal of flaps and splitting bone.

SURGICAL CHISELS



2001  
13K/TG 3.5/3mm



3.5mm

2001T  
13K/TG TITANIUM

3mm



2001  
13K/TG 3.5/4.5mm



3.5mm

2002T  
13K L/R TITANIUM

4.5mm



4mm

2002L  
13KL 4/6mm  
TITANIUM AVAILABLE

6mm

KRAMER NEVINS SINUS LIFT



2003  
KN1/2



2004  
C36/37



2005  
44/45



For scraping out and reshaping bone.  
For removal of flaps and splitting bone.

KRAMER NEVINS SINUS LIFT



2006  
TG



2007  
SOLT 3/4



2008  
TG-O



2010 OCH-1



2011 OCH-2

BONE CHISELS



2009  
29K

2019  
KN-3

2018  
KN-55

2012  
FEDI-4



2013  
OCH-3



2014  
OCH-4



2015  
FEDI-1



2016  
FEDI-2



2017  
FEDI-3



2020



For mixing materials and general laboratory use.



# Implant Instruments D

## IMPLANT KIT 27 INS (Recommended)

1. 1105 1EZ COLOR PROBE CP15UNC
2. 1923T 6EZ PERIOTOME PT2 TITANIUM WITH SERRATION
3. 1924T 6EZ PERIOTOME PT1 TITANIUM WITH SERRATION
4. 2033N OCT JOVANOVIC 3MM BONE CONDENSER
5. 2034S 6EZ BONE PACKER IBP3/5 SERRATION
6. 2071 FORCEP KERISON
7. 2094 BONE MIXING BOWL
8. 3506 6EZ PERIOSTEALS MOLT9
9. 3751 ELEVATORS LUXATION 2MM STD
10. 4103 6EZ SURGICAL CURETTES LUCAS86
11. 4111 6EZ SURGICAL CURETTES MOLT2/4
12. 4618 STAINLESS RULER
13. 4625 BONE CALIPER DOUBLE SCALE
14. 5505TN MICRO NEEDLE HOLDER T/C
15. 5917 GUM SCISSORS LA GRANGE 11.5CM
16. 5920 GUM SCISSORS IRIS CURVED 11.5CM
17. 6420 RETRACTORS UNIVERSITY OF MINNESOTA
18. 6428 RETRACTORS BISHOP
19. 6502 FORCEP ADSON 1X2 12CM
20. 6717 SUCTION TUBES 4MM
21. 6905 SCALPEL HANDLE 5
22. 6923 SCALPEL HANDLE ADJ
23. 7018 SYRINGES ASPIRATED
24. 7425 BONE TRAY
25. 7505 CASSETTE 20 HOLD
26. 7505-1 CASSETTE CLIP SCISSOR
27. 7505-2 CASSETTE CLIP SUCTION



Black Titanium  
Knives  
Surgical Chisels  
Osteotomes  
Sinus Lift Instruments  
Osteotomes  
Implant Instruments

Implant Scalers  
Tissue Pliers  
Castroviejo  
Accu Luxating  
Titanium  
Implant Scalers

HANDLE OPTION

\* 6EZ D/E OR S/E

## BLACK TITANIUM INSTRUMENTS

Our black titanium instruments are enhanced in strength and durability up to 3 times more than the ordinary tips. It is highly hygienic as it protects the instruments from virus infection. It is easy to clean and requires less sharpening than other instruments.

The titanium provides the surface of the instrument to be a special non-stick design.



1601T  
GR1/2  
BLACK TITANIUM

## BLACK TITANIUMS



1201T  
H6/7  
BLACK TITANIUM



1203T  
H5/33  
BLACK TITANIUM



1208T  
204S  
BLACK TITANIUM



1226T  
NV-2  
BLACK TITANIUM



1322T  
MC13/14S  
BLACK TITANIUM



1415T  
COL4R/4L  
BLACK TITANIUM



1416T  
COL13/14  
BLACK TITANIUM



1422T  
BH5/6  
BLACK TITANIUM



1601T  
GR1/2  
BLACK TITANIUM



1606T  
GR11/12  
BLACK TITANIUM



1607T  
GR13/14  
BLACK TITANIUM



1608T  
GR15/16  
BLACK TITANIUM



1609T  
GR17/18  
BLACK TITANIUM



1610T  
AF1/2  
BLACK TITANIUM



1614T  
AF11/12  
BLACK TITANIUM



1615T  
AF13/14  
BLACK TITANIUM



1620T  
SUB1/2  
BLACK TITANIUM



1621T  
SUB11/12  
BLACK TITANIUM



1622T  
SUB3/4  
BLACK TITANIUM





1812  
KOR1/2



1812M  
KOR1/2M



1812L  
KOR1/2L



1813  
KCK-3



1816D  
KM3/4



1901  
PR1/2



1902  
KRA-1



1903  
KRA-2



1904



SURGICAL CHISELS



2001  
13K/TG 3.5/3mm



3.5mm

2001T  
13K/TG TITANIUM

3mm



2002  
13K/TG 3.5/4.5mm



3.5mm

2002T  
13K L/R TITANIUM

4.5mm



4mm

2002L  
13KL 4/6mm  
TITANIUM AVAILABLE

6mm



Stainless OCH-1

2010 OCH-1





STAINLESS OCH-2

2011 OCH-2



For bone spreading for implant.



**1930**  
2mm STR



**1930C**  
2mm CVD



## SINUS LIFT INSTRUMENTS

Allows for transport of the bone graft chips from implant cavity or for condensing them apically.  
For separation of sinus lining from bone.



Commonly used implant instruments.



2029 CS1ABANCA BONE CURETTE 7mm



2030 BSJO 5 JOVANOVIC 5mm BONE SCOOP



2031 P 25 PERIOSTEAL



2032 LABANCA BONE PLUGGERS 3X1.5mm



2033 JOVANOVIC BONE PLUGGERS 3/3mm



2033A  
BONE GRAFT 1.0/1.7mm TITANIUM  
8-11-13-15-18



2033B  
BONE GRAFT 1.7/2.0mm TITANIUM  
8-11-13-15-18



2033C  
BONE GRAFT 2.4/2.7mm TITANIUM  
8-11-13-15-18



2033N  
BONE GRAFT 3.3/4mm  
8-11-13-15-18

2033NS  
BONE GRAFT 3.3/4mm  
SERRATED  
8-11-13-15-18



2033NT  
BONE GRAFT 3.3/4mm  
TITANIUM  
8-11-13-15-18

2033NTS  
BONE GRAFT 3.3/4mm  
TITANIUM SERRATED  
8-11-13-15-18



2034S BONE PACKER IBP 3/5mm SERRATED



BONE PACKER  
SERRATED



Commonly used implant instruments.



1923T  
PT-2 TITANIUM TG  
A/V Regular



1924T  
PT-1 TITANIUM TG  
A/V Regular



1925T  
PT-3 TITANIUM Angel S/E



1926T  
PT-4 TITANIUM STR S/E



1927T  
PT-5 TITANIUM TG CVD S/E



1928  
MS-1 PERIOTOME



2043T  
PPAEL TITANIUM  
REGULAR AVAILABLE



2044T  
PPAELA TITANIUM  
REGULAR AVAILABLE




2035T BONE CARRIER TITANIUM  
2035 BONE CARRIER AVAILABLE



2038T BONE GRAFT/CARRIER TITANIUM 2-4-6-8-10

OSTEOTOMES - Used in surgical procedures for preparation of implants.



**2050S**  
OSTM 3.2mm STR  
8-10-13-15-18mm



**2051S**  
OSTM 4.0mm STR  
8-10-13-15-18mm



**2052S**  
OSTM 1.8-3.3mm STR  
8-10-13-15-18mm




**2053S**  
OSTM 2.5-4.0mm STR  
8-10-13-15-18mm



**2054S**  
OSTM 3.0-4.7mm STR  
8-10-13-15-18mm



**2055S**  
OSTM 4.0mm STR  
8-10-13-15-18mm



**2056S**  
OSTM R0.38-1.8mm STR  
8-10-13-15-18mm



STRAIGHTH

**2050S**  
OSTM 3.2mm STR  
8-10-13-15-18mm

Osteotome - Surgical instruments used to cut or divide bone.



**2057**  
OSTM 4mm STR  
7-10-13-15-18mm



**2058**  
OSTM 4.0mm CVD  
7-10-13-15-18mm




**2059**  
OSTM 7.5mm STR  
7-10-13-15-18mm



**2060**  
OSTM 7.5mm CVD  
7-10-13-15-18mm



**2061**  
OSTM 6.0mm  
JOVANOVIC



**2064**  
SURGICAL CHISEL  
3mm 4EZ S/CUT



**2064D**  
SURGICAL CHISEL  
3mm 4EZ D/CUT



**2066**  
SURGICAL CHISEL  
9mm 4EZ S/CUT



**2065**  
SURGICAL CHISEL  
6mm 4EZ S/CUT



**2066D**  
SURGICAL CHISEL  
9mm 4EZ D/CUT



**2065D**  
SURGICAL CHISEL  
3mm 4EZ D/CUT



**2067**  
FORCEP BONE GRAFT  
STR 19cm



**2068**  
FORCEP BONE GRAFT  
CVD 19cm



**2067**  
FORCEP BONE GRAFT  
STR 19cm

\* STR : Straight      CVD : Curved



Commonly used implant instruments.



2071 FORCEP KERISON 3mm x 75mm  
2072 FORCEP KERISON 4mm x 75mm



2073 ENDO FORCEP 75mm  
2074 ENDO FORCEP 95mm  
2075 ENDO FORCEP 110mm

2079 BONE SYRINGE 3mm  
2080 BONE SYRINGE 3.5mm  
2081 BONE SYRINGE 4mm  
2082 BONE SYRINGE 7mm



2083  
BONE MORSELIZER




2087 BONE CRUSHER



2091 BONE COLLECTOR



2093 BONE COLLECTOR CHISEL STR  
2093A BONE COLLECTOR CHISEL CVD



2093B BONE COLLECTOR STR BLADE ONLY  
2093AB BONE COLLECTOR CVD BLADE ONLY



2092 BONE COLLECTOR WITH FILTER

2092F BONE COLLECTOR FILTER ONLY



2094  
BONE MIXING BOWL



Commonly used implant instruments.



4613  
CASTROVIEJO  
STR 8.5cm



4613L  
CASTROVIEJO  
STR 17cm



4613LA  
CASTROVIEJO  
CVD 17cm



4613RA  
CASROVIEJO  
17cm



4614  
Ridge Mapping Bone width 0 ~25mm  
14cm



4618  
RULLER  
15cm / 6in



4625  
CALIPER DOUBLE  
15cm



6420A  
LINDO  
SERRATED



6435  
BAYONET 15cm



Commonly used implant instruments.



6931 PARA ANTERIOR 5mm



6932 PARA POSTERIOR 5mm



4411  
MICRO FRIEDMAN 15  
14.5cm



4412  
MICRO FRIEDMAN 90  
14.5cm



4413  
MICRO FRIEDMAN 45  
14.5cm



6436  
WEIDER TONGUE  
30mm

6436L  
WEIDER TONGUE LARGE  
35mm

CHALAZION MUCOCELE



6442  
10 X 10 mm



6443  
11 X 16 mm



6444  
12 X 21 mm



6445  
LIP RETRACTOR



7425  
TRAY IMPLANT



7426  
TRAY IMPLANT

The ergonomic hold - The design of tissue and cotton pliers prevents slipping while reducing the weight over 20%.

### UNIQUE GRIP DESIGN



6501  
ADSON 12cm



6501L  
ADSON 15cm



6501A  
ADSON 12cm



6502  
ADSON 12cm



6502L  
ADSON 15cm



6502B  
ADSON 12cm

## Commonly used implant instruments.



6511  
SUTURE 16cm



6511N  
SUTURE 16cm



6512  
GOLDMAN-FOX  
(CRANE-KAPLANE)  
15cm



6513  
GOLDMAN-FOX  
(CRANE-KAPLANE)  
15cm



6522  
SELF-LOCKING 16cm



6525  
TISSUE FORCEPS T/C  
GERALD 18cm




6526  
TOMOR 16cm



6526L  
TOMOR 20cm







**5309**  
BABY-CRILE (PEAN)  
STR 14cm



**5926 CVD**  
CASTROVIEJO CVD  
10cm/4"

**5926 STR**  
CASTROVIEJO STR  
10cm/4"



**5929**  
BARRAQUER  
14cm

AVALBLE ROUND HANDLE  
**5929C** CURVED



**5944**  
CASTROVIEJO SCISSORS  
45° T/C



**4223M**  
130g  
MALLET



**4223**  
235g  
MALLET



**7018**  
SYRINGE 1.8ml



**6923**  
SCALPEL 180°  
16.2cm

For holding and guiding the needle securely during suturing.



**5513**  
CASTROVIEJO T/C  
16cm 5 1/2"



**5514**  
CASTROVIEJO T/C  
18cm 7"



**5515**  
CASTROVIEJO T/C  
CURVED 16cm 5 1/2"



**5516**  
CASTROVIEJO T/C  
CURVED 18cm 7"



**5409C**  
CASTROVIEJO T/C  
16CM CVD



**5409**  
CASTROVIEJO  
16cm STR



**5525**  
STRAIGHT WITH SPOON  
15cm



**5526**  
CURVED WITH SPOON  
15cm



**5527**  
CURVED  
14cm



**5531**  
CASTROVIEJO STRAIGHT T/C  
CVD 45°



**5505TN**  
MICRO NEEDLE HOLDER T/C  
15cm

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.

BLUE TITANIUM WITH SERRATED



3769 ACCU LUX STRAIGHT



3770 ACCU LUX CURVED



3771 ACCU LUX SPADE



3772 ACCU LUX DISTAL



3773 ACCU LUX MESIAL



3774 ACCU LUX-LARGE DISTAL



3775 ACCU LUX-LARGE MESIAL



TITANIUM INSTRUMENTS



3506T  
TITANIUM MOLT-9



3509T  
TITANIUM PR-3



3514T  
TITANIUM 24G



3516T  
TITANIUM BUSER



3521T  
TITANIUM HF-20




3523  
SELDEN-23



3526  
WOODSON-1 MODIFY



# Operative E



## COMPOSITE SET UP 19 INS (Recommended)

- |     |       |     |                             |
|-----|-------|-----|-----------------------------|
| 1.  | 1006  | 3EZ | EXPLORERS EXP 5             |
| 2.  | 2231  | 4EZ | EXCAVATORS 32L              |
| 3.  | 2501T | 4EZ | CARVERS IPC, TITANIUM       |
| 4.  | 2702  | 4EZ | CARVING KNIVES 21           |
| 5.  | 2814  | 4EZ | PICH                        |
| 6.  | 2901T | 6EZ | PLASTIC FILLING(THIN) 1, TN |
| 7.  | 2904T | 6EZ | PLASTIC FILLING(THIN) 4, TN |
| 8.  | 2910T | 6EZ | PLASTIC FILLING TUFTS4, TN  |
| 9.  | 3103  | 6EZ | GINGIVAL RETRACTORS GRM2    |
| 10. | 3108  | 6EZ | CORD PACKER 113             |
| 11. | 4310  |     | AMALGAM CARRIERS M/J        |
| 12. | 4706  |     | FORCEP PAPER                |
| 13. | 6517  |     | FORCEP LONDON-COLLEGE       |
| 14. | 6911  | 6EZ | MIRROR HANDLE #6 C/S        |
| 15. | 6951  |     | #5 MIRROR C/S FRONT SURFACE |
| 16. | 7213  |     | #207 RUBBER DAM CLAMPS      |
| 17. | 7217  |     | #27 RUBBER DAM CLAMPS       |
| 14. | 7233  |     | #210 RUBBER DAM CLAMPS      |
| 15. | 7236  |     | IVORY LIGHT WEIGHT 17CM     |
| 16. | 7502  |     | CASSETTES 8-16 HOLD         |
| 17. | N     |     | #1 RUBBER DAM CLAMP         |
| 18. | N     |     | #14 RUBBER DAM CLAMP        |
| 19. | N     |     | #9 RUBBER DAM CLAMP         |

## HANDLE OPTION

- \* 4EZ D/E OR S/E
- \* 6EZ D/E OR S/E

Chisels  
 Excavators  
 Plugger/Condensers  
 Carvers  
 Carving Knives  
 Burnishers  
 Gingival Retractors  
 Cord Packers  
 PKT Waxing  
 Amalgam Carriers

- \* Hatchets : For cavity preparation and removing hard caries.
- \* Hoes : For forming line angles on anterior area.
- \* Chisels : For refining the cavity preparation.
- \* Margin Trimmer : For producing proper bevel on enamel margins.

HATCHETS



CHISELS



MARGIN TRIMMERS




HOE






For removing debris and decay from cavities.



For removing debris and decay from cavities.




For packing and filling materials into cavity preparations.



\* T A/V : Available Titanium

For packing and filling materials into cavity preparations.






For carving anatomical features, smoothing surfaces and trimming excess materials.




ANDREWS




CLAPPISON



HOLLENBACK CLEOID




HOLLENBACK




SHOSHAN



THARP



For trimming excess filling materials.



For condensing and smoothing carved amalgam.






■ M: Mini ■ B: Boghosian ■ Titanium Available




For composite placing and shaping.

THIN COMPOSITE



UNIVERSAL




- M: Mini ■ B: Boghosian ■ T/G: Titanium-Gold Available
- B: BEQ(Quetin) ■ T: Tufts ■ WDS: Woodson ■ Titanium-Available

For composite placing and shaping.

ANTERIOR



POSTERIOR



- \* Gingival Retractors : For protecting tissue during cavity preparations and placement of composite fillings.
- \* Cord Packer : For protecting tissue of gingival area by putting the thread between teeth and gingival area.

GINGIVAL RETRACTORS



CORD PACKERS



For shaping, carving and flowing of molten wax.



3203  
3R



3201  
1



3202  
2



3204  
3



3205  
4



3206  
5



3211  
GOLD FOGEL-1



For carrying and dispensing amalgam materials.



4301  
SMALL 1.5mm



4302  
MEDIUM 2mm



4303  
LARGE 2.5mm



4304  
JUMBO 3mm



4305 4306 4307

REMOVABLE TIP



4311  
AMALGAM WELL SMALL



4308  
-0.15mm- SMALL/MEDIUM -2mm-



4310  
-2mm- MEDIUM/JUMBO -2.5mm-



4310 EZ  
-2mm- EZ MEDIUM/JUMBO -2.5mm-



2404T  
PLG-3 TITANIUM



2407T  
PLGW-3 TITANIUM



2413T  
PLGO-1 TITANIUM



2415T  
PLGO-3 TITANIUM



2415TR  
PLGO-TR TITANIUM



2420T  
PLG S/M TITANIUM



2421T  
PLG ML TITANIUM



2422T  
KTO TITANIUM



2501T  
IPC TITANIUM



2502T  
IPCOA TITANIUM



2801T  
B-18 TITANIUM



2802T  
B26/27S TITANIUM



2803T  
B27/29 TITANIUM



2806T  
B21 TITANIUM



2807T  
B21L TITANIUM



2820T  
B28/29 TITANIUM



2821T  
B29/34 TITANIUM



2825T  
RU-1 TITANIUM



2826T  
RU-2 TITANIUM



2827T  
RU-3 TITANIUM



2901T  
PFIT-1 TITANIUM



2902T  
PFIT-2 TITANIUM



2903T  
PFIT-3 TITANIUM



2904T  
PFIT-4 TITANIUM





2910T  
PFI-4 TITANIUM



2920T  
WDS-2 TITANIUM



2921T  
PFI CM-11 TITANIUM



2922T  
PFIM-12 TITANIUM



2923T  
PFIT-5 TITANIUM



3001T  
PF18A TITANIUM




3002T  
PFIT-11 TITANIUM

# Endodontic F

ENDODONTIC SET UP 16 INS (Recommended)

1. 1006 3EZ EXPLORERS EXD 5
2. 1026 3EZ EXPLORERS DG16
3. 1104 3EZ EXPLORERS CP12
4. 2230 4EZ EXCAVATORS 31L
5. 3308 3EZ PLUGGERS/CONDENSERS 5/7
6. 3410 3EZ ENDODONTIC SPREADERS D11T
7. 3416 4EZ ENDODONTIC HEAT CARRIERS MA57
8. 6520 FORCEP SELF-LOCKING 16CM
9. 6911 6EZ MIRROR HANDLE #6 C/S
10. 6951 #5 MIRROR C/S FRONT SURFACE
11. 7236 IVORY LIGHT WEIGHT 17CM
12. 7502 CASSETTES 16 HOLD DOUBLE DECK
13. N #1 RUBBER DAM CLAMP
14. N #2 RUBBER DAM CLAMP
15. N #9 RUBBER DAM CLAMP
16. N #26 RUBBER DAM CLAMP




Pluggers & Condensers  
Spreaders  
Pliers  
Rubber Dam Clamps  
& Accessories

HANDLE OPTION

\* 3EZ D/E OR S/E

\* 4EZ D/E OR S/E


For vertical pressing and compacting of filling materials.



- \* For filling solid medicine materials into the endodontic area during lateral condensation.
- \* The function of spreaders is similar to Pluggers and Condensers.
- \* For lateral pressing and compacting of filling materials.



RCS SET



HEAT CARRIERS







4627  
ENDO FINGEL ROLLER



6515  
SELF-LOCKING 15cm



6517L  
LONDON-COLLEGE  
(SELF-LOCKING)  
15cm



6523  
SELF-LOCKING

Rubber Dam Clamp – Stabilizes and holds the dam secure.

ANTERIOR



PREMOLARS



LOWER MOLARS



Rubber Dam Clamp – Stabilizes and holds the dam secure.

UPPER MOLARS



AVAILABLE : 31A, 51A

MOLARS - UNIVERSAL



AVAILABLE : 51, 138, 26A, 28A,  
30A, 203, 204



**7236**  
 IVORY  
 LIGHT WEIGHT  
 17cm



**7237**  
 BREWER  
 17.5cm



**7238**  
 STOKE  
 17.5cm



**7239**  
 AINSWORTH  
 17cm



**7240**  
 IVORY  
 17cm



**7241S** RUBBER DAM FRAME 9cm  
**7241M** RUBBER DAM FRAME 13cm  
**7241L** RUBBER DAM FRAME 15cm

S:9cm



# Surgical G

## SURGICAL SET UP 13 INS (Recommended)

1. 1812 6EZ SURGERY KNIVES KOR1/2
2. 2001 6EZ SURGICAL CHISELS 13K/TG
3. 3509 6EZ PERIOSTEALS PR3
4. 3513 6EZ PERIOSTEALS ALLEN9
5. 3514 6EZ PERIOSTEALS 24G
6. 4107 6EZ SURGICAL CURETTES LUCAS87
7. 4411 RONGEURS MICRO FRIEDMAN 15' (14.5CM)
8. 4618 STAINLESS STEEL RULER
9. 5505TN T/C MICRO CRILE-WOOD STR 15CM
10. 5905 KELLEY SCI, 16CM
11. 6420 RETRACTORS MINNESOTA
12. 6501 PLIERS ADSON 12CM
13. 6905 SCALPEL HANDLE 5



Periosteals  
Surgical Currettes(Spoon)  
Bone File/Chisels  
Bone Rongeurs  
Root Splinters  
Haemostatic Forceps  
Needle Holders  
Gum Scissors  
Tungsten Carbide Inserts

Raspatories  
Skin Hooks  
Retractors  
Scalpel Handles  
Syringes

### HANDLE OPTION

\* 6EZ D/E OR S/E

\* Elevator Small OR Elevator Large

For retracting periosteum after removing gingival tissue.



3501  
BENNETT



3502  
P-149



3503  
HOPKINS-1



3504  
HOPKINS-2



3505  
MEAD-3



3505M  
MEAD-3MINI



3506  
MOLT-9



3506T  
TITANIUM MOLT-9



3507  
GF-14

For retracting periosteum after removing gingival tissue.





For retracting periosteum after removing gingival tissue.



3514  
24G



3514T  
TITANIUM 24G



3515  
BENQUE



3516  
BUSER



3516T  
TITANIUM BUSER



3517  
KN-1



3518  
KN-152



3520  
HOURIGAN-2



For retracting periosteum after removing gingival tissue.



For retracting periosteum after removing gingival tissue.



3531  
P25



3532  
P25L



3533  
P25XL

For curettage and cyst removal.



4101  
84 LUCAS



4102  
85 LUCAS



4103  
86 LUCAS



4107  
87 LUCAS



4108  
88 LUCAS

SURGICAL CURETTES (SPOON)

For curettage and cyst removal.



4104  
8 MILLER



4105  
9 MILLER



4106  
10 MILLER



4109  
11 MILLER



4110  
12 MILLER



4111  
MOLT 2/4



4112  
MOLT-2



4113  
MOLT-4



4114  
MOLT-5L



4115  
MOLT-5R



- \* Bone Chisels : For final smoothing of bone by pushing out.
- \* Bone Files : For splitting and removing bone by pulling out.



4219  
BF1X

SINGLE BEVEL



4201 1 GARDNER    4202 2 GARDNER    4203 5 GARDNER    4204 5R GARDNER    4205 6 GARDNER    4206 9 GARDNER    4207 1 CHANDLER    4208 2 CHANDLER    4209 4 CHANDLER

BI BEVEL



4210 1 WAKEFIELD    4211 2 WAKEFIELD    4212 18 MCFARLAND    4213 7 CHANDLER    4214 8 CHANDLER    4215 5 CHANDLER

4216  
BF45  
MILLER

4217  
BF1W  
WAHL

4218  
BF2W  
WAHL



4219  
BF1X  
MILLER

4220  
BF2X  
MILLER

4221  
BF3X  
MILLER

4222  
BF21  
MILLER



4223M  
130g  
MALLET



4223  
235g  
MALLET



4224  
MEAD-2  
MALLET

Gouges : For scooping the bone.



4501  
PARTSCH  
17cm



4502  
FREER  
16cm



4503  
FREER  
16cm



4504  
LUCAS  
16cm



4505  
LUCAS  
16cm



4506  
LAMBOTTE  
12.5cm



4507  
14cm



4508  
PARTSCH  
14cm



4509  
PARTSCH  
17cm



4510  
PARTSCH  
14cm

For tooth extraction as a whole by trimming and recontouring alveolar bone.  
For gross tissue removal.



**4401**  
BLUMENTHAL  
15cm



**4402**  
BLUMENTHAL  
45° 14.5cm



**4403**  
BLUMENTHAL  
90° 14.5cm



**4404**  
MEAD 1A  
14cm



**4405**  
BEYER 18cm



**4406**  
BEYER MINI 15.5cm



**4407**  
LUER 15cm



**4408**  
FRIEDMAN 14cm



**4409**  
MINI-FRIEDMAN 15  
12cm



**2071** FORCEP KERISON 3mm  
**2072** FORCEP KERISON 4mm



**4411**  
MICRO FRIEDMAN 15  
14.5cm



**4412**  
MICRO FRIEDMAN 90  
14.5cm



**4413**  
MICRO FRIEDMAN 45  
14.5cm

ROOT SPLINTER

For extracting root tips of the tooth where normal forceps cannot reach.



5105  
ROOT FRAGMENTS  
STRAIGHT



5106  
45° Degree



5107  
Available  
90° Degree



PEET SHARP



5108  
15° Degree



5109  
45° Degree



5110  
90° Degree





**5201**  
WITHOUT RATCHET  
14.5cm



**5202**  
WITHOUT RATCHET  
14.5cm



**5203**  
CURVED  
17.5cm



**5204**  
CURVED  
21cm



**5205**  
WITHOUT RATCHET  
18cm



**5206**  
WITHOUT RATCHET  
18cm



**5207**  
WITHOUT RATCHET  
20cm



**5208**  
WITHOUT RATCHET  
20cm

Haemostatic Forceps – Locking forceps for seizing the end of a blood vessel to control hemorrhage.



5301C  
CURVED

5301  
BABY-MOSQUITO  
STR 10cm



5302C  
CURVED

5302  
MICRO-HALSTED  
STR 12.5cm



5303C  
CURVED

5303  
HALSTED-MOSQUITO  
STR 12.5cm



5305C  
CURVED

5305  
MICRO-HALSTED  
STR 12.5cm



5307  
PEAN 12.5cm



5309C  
CURVED

5309  
BABY-CRILE (PEAN)  
STR 14cm



5310C  
CURVED

5310  
CRILE  
STR 14cm



5311C  
CURVED

5311  
BABY-CRILE (KOCHER)  
STR 14cm



53012C  
CURVED

5312  
CHILE (KOCHER)  
16cm



53013C  
CURVED

5313  
HALSTED  
STR 18.5cm



53015C  
CURVED

5315  
HALSTED  
STR 18.5cm



5316  
BABY-MIXTER  
CURVED 18cm

For holding and guiding the needle securely during suturing.



For holding and guiding the needle securely during suturing.





High-tech welding of tungsten carbide greatly improves the durability and serration. Exceptionally upgrades the performance.



**5504**  
OLSON-HEGAR  
15cm



**5505** CRILE-WOOD 15cm  
**5505L** CRILE-WOOD 18cm



**5505TN**  
MICRO  
15cm



**5506**  
CRILE-WOOD CURVED  
15cm



**5507** RYDER 15cm  
**5508** RYDER 18cm



**5509**  
DERF  
11.5cm



**5510** HEGAR-BAUMGARTEN 14cm  
**5511** HEGAR-BAUMGARTEN 13cm  
**5512** HEGAR-BAUMGARTEN 13cm



**5517**  
ARRUGA  
STR 12cm

**5518**  
ARRUGA  
CURVED  
12cm



**5519**  
WEBSTER  
11.5cm



**5521**  
MAYO-HEGAR  
STRAIGHT  
18cm

**5523**  
MAYO-HEGAR  
CURVED  
16cm

**5521L**  
MAYO-HEGAR  
STRAIGHT  
20cm

**5524**  
MAYO-HEGAR  
CURVED  
18cm

**5520**  
MAYO-HEGAR  
STRAIGHT  
16cm

For holding and guiding the needle securely during suturing.



**5513**  
CASTROVIEJO T/C  
16cm 5 1/2"



**5514**  
CASTROVIEJO T/C  
18cm 7"



**5515**  
CASTROVIEJO T/C  
CURVED 16cm 5 1/2"



**5516**  
CASTROVIEJO T/C  
CURVED 18cm 7"



**5409C**  
CASTROVIEJO T/C  
16CM CVD



**5409**  
CASTROVIEJO  
16cm STR



**5525**  
STRAIGHT WITH SPOON  
15cm



**5526**  
CURVED WITH SPOON  
15cm



**5527**  
CURVED  
14cm



**5531**  
CASTROVIEJO STRAIGHT T/C  
CVD 45°



**5505TN**  
MICRO NEEDLE HOLDER T/C  
15cm

For cutting threads and gum.



For cutting threads and gum.



**5911**  
GOLDMAN-FOX  
12.5cm



**5912**  
MAYO  
14.5cm

**5913C**  
MAYO CVD  
14.5cm



**5916C**  
GOLDMAN-FOX CVD  
12.5cm



**5916TC**  
GOLDMAN-FOX CVD  
12.5cm



**5917C**  
LA GRANGE CVD  
11.5cm



**5917TC**  
LA GRANGE CVD  
11.5cm



**5919**  
IRIS  
11.5cm



**5920C**  
IRIS CVD  
11.5cm



For cutting threads and gum.



5923  
SUTURE  
12cm/4 3/4"



5925  
CROWN  
12cm/4 3/4"  
SAME AS #4602



5926 CVD  
CASTROVIEJO CVD  
10cm/4"



5928  
GENERAL  
14.5cm/5 3/4"

5926 STR  
CASTROVIEJO STR  
10cm/4"



5929  
BARRAQUER  
14cm



5930  
GENERAL  
14.5cm/5 3/4"



5931  
BANDAGE  
14.5cm/5 3/4"



5935  
16cm



5940  
CASTRO CVD



5941  
CASTRO CVD



5944  
CASTROVIEJO SCISSORS  
45° T/C

AVAILABLE ROUND HANDLE  
5929C CURVED

For cutting threads and gum.



6001  
METZENBAUM  
14.5cm



6002C  
METZENBAUM CVD  
14.5cm



6003  
METZENBAUM  
14.5cm



6004  
METZENBAUM CVD  
14.5cm



6005  
METZENBAUM  
18cm



6006  
GOLDMAN-FOX  
13cm



6007C  
GOLDMAN-FOX CVD  
13cm



6008  
WIRE CUTTING  
12cm



6009  
IRIS STR  
11.5cm



6010  
IRIS CVD  
11.5cm



6011  
DEAN  
16.5cm



6012  
LOCKLIN  
16.5cm



6013  
LOCKLIN  
16.5cm



6014  
KELLEY  
16cm

\* Skin Hooks : For hooking skin during surgery.

\* Retractors : For holding back the lips and cheeks.



Skin Hooks : For hooking skin during surgery.

Retractors : For holding back the lips, cheeks, and tongue.



6301  
FRAZIER  
13cm



6305  
JOSEPH  
16.5cm



6312  
GILLIES  
18cm



6325  
SHARP



6327  
SHARP



6332  
16cm



6333



6334



6335  
16cm



6336  
STANDELL-STILLE  
18cm



6337  
BLUNT



6338  
TOOTHED



6339  
LANGENBECK-GREEN  
16cm



6340



6341



6346  
TR28K-L



6346  
TR28K-R



25X10mm  
22X10mm

6342  
FARABEUF 12.5cm

28X12mm  
32X12mm



23X16mm  
20X16mm

6343  
FARABEUF 15cm

24X16mm  
28X16mm



29mm  
26mm  
22mm

6344  
ROUX  
16.5cm

28mm  
38mm  
45mm



Retractors : For holding back the lips, cheeks, and tongue.

Mouth Gag : To prop open the mouth and prevent it from closing.



Retractors : For holding back the lips, cheeks, and tongue.



6432  
RETRACTORS



6433  
RETRACTORS LARGE



6434  
RETRACTORS SMALL



6436  
WEIDER TONGUE  
30mm



6436L  
WEIDER TONGUE LARGE  
35mm



6435  
BAYONET 15cm



6437  
SPANDEX LARGE 12.5cm



6438  
SPANDEX SMALL 10cm



6439  
RETRACTOR D/E  
14cm



6440  
RETRACTOR S/E  
9cm



6442  
10 X 10 mm



6443  
11 X 16 mm




6444  
12 X 21 mm



6445  
LIP RETRACTOR

CHALAZION MUCOCELE

Scalpel handle : For holding a mess blade.



SYRINGES – An instrument for injecting anesthetics from a cartridge.

3-WAY SYRINGES – Supplies air, water, or a mist through separate interval channels created by combining pressurized air with water flow.



7001  
INTRALIGMENTARY  
SYRINGE 1.8ml



7015  
SYRINGE 1.8ml



7015L  
SYRINGE 1.8ml



7003  
SYRINGE 1.8ml



7017  
PAROJECT



7018  
SYRINGE 1.8ml



7012  
SYRINGE 1.8ml



7018L  
SYRINGE 2.2ml



7013  
SYRINGE 1.8ml



7020  
PAROJECT



7019  
WATER SYRINGE 3WAY





# Elevators H

## SURGICAL SET UP 13 INS (Recommended)

1. 1812 6EZ SURGERY KNIVES KOR1/2
2. 2001 6EZ SURGICAL CHISELS 13K/TG
3. 3509 6EZ PERIOSTEALS PR3
4. 3513 6EZ PERIOSTEALS ALLEN9
5. 3514 6EZ PERIOSTEALS 24G
6. 4107 6EZ SURGICAL CURETTES LUCAS87
7. 4411 RONGEURS MICRO FRIEDMAN 15' (14.5CM)
8. 4618 STAINLESS STEEL RULER
9. 5505TN T/C MICRO CRILE-WOOD STR 15CM
10. 5905 KELLEY SCI, 16CM
11. 6420 RETRACTORS MINNESOTA
12. 6501 PLIERS ADSON 12CM
13. 6905 SCALPEL HANDLE 5



### HANDLE OPTION

- \* Elevator Small
- \* Elevator Large

Elevators  
Root Tip Picks/Teasers  
Accu Luxating

ELEVATORS

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.

All Handles are Unique Design  
Shown actual size.



**3701**  
HEID BRINK 1



**3702**  
HEID BRINK 2



**3703**  
SELDIN 4S



**3704**  
SELDIN 12M



**3705**  
SELDIN 34



**3706**  
SELDIN 34S



**3707**  
SELDIN H40



**3708**  
HEID BRINK 41



**3709**  
APEXO 80



**3710**  
APEXO 81



**3711**  
APEXO 81E



**3712**  
APEXO 301



**3713**  
APEXO 302



**3714**  
SELDIN 2S



**3715**  
APEXO 303



**3716**  
CRANE PICK 8



**3717**  
SELDIN 1L



**3718**  
SELDIN 1R



**3719**  
SELDIN 4L



**3720**  
SELDIN 4R



**3721**  
CRYER M 2S



**3722**  
CRYER M 26

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.

All Handles are Unique Design  
Shown actual size.



3723  
CRYER L 31



3724  
CRYER L 32



3725  
CRYER S 30A



3726  
CRYER S 31A



3727  
CRYER M 44



3728  
CRYER M 45



3729  
CRYER L 46L



3730  
CRYER L 46R



3731  
BARRY 50



3732  
BARRY 51



3734  
MILLER APEXO  
71



3735  
MILLER APEXO  
72



3736  
MILLER APEXO  
73



3737  
MILLER APEXO  
74



3738  
MILLER APEXO  
75



3739  
APICAL  
302 3L



3740  
APICAL  
301 2S



3741  
APICAL  
303 4R



3742  
FLOHR L  
SHARP END



3743  
FLOHR STR  
SHARP END



3744  
FLOHR R  
SHARP END

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



3748  
BURNARD  
2 CVD



3747  
BURNARD  
3L



3748  
BURNARD  
2 STR



3749  
BURNARD  
2 OUT



3750  
BURNARD  
4R



3751  
LUXATION  
2mm STRAIGHT



3752  
LUXATION  
3mm STRAIGHT



3753  
LUXATION  
5mm STRAIGHT



3754  
LUXATION  
2mm CURVED



3755  
LUXATION  
3mm CURVED



3756  
LUXATION  
5mm CURVED



3757  
CRYER UK  
30L



3758  
CRYER UK  
30R



3759  
COUPLAND 3  
4.5mm



3760  
COUPLAND 2  
3.5mm



3761  
COUPLAND 1  
2.5mm



3762  
HYLIN 3  
4mm



3763  
HYLIN 1  
2mm



3764  
HYLIN 2  
3mm



3776  
LUX-TOME 3mm



3777  
LUX-TOME 4mm



3778  
LUX-TOME 5mm



For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



3765  
LINDO 3mm



3766  
LINDO 4mm



3767  
E 46



3768  
E 46R



3830  
LINDO LEVIN-E93  
SERRATED



3831  
LINDO LEVIN-E92  
SERRATED



3841  
COGSWELL B



3846  
LUXATION  
3mm CURVED



3842  
LARGE WING 1



3843  
LARGE WING 2



3844  
LARGE WING 3



3845  
LARGE WING 4

SMALL HANDLE AVAILABLE



3852  
POTTS 6



3853  
POTTS 7



4606  
CROWN SPREADER

ELEVATORS

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



3747S  
BURNARD  
3 S



3748S  
BURNARD  
2S STR



3750S  
BURNARD  
4S



3803  
FLOHR 3R



3802  
FLOHR 1STR



3801  
FLOHR 2L



3804  
KOPP  
2L



3805  
KOPP  
1STR



3806  
KOPP  
3R



3810  
FLOHR  
4mm



3811  
FLOHR  
5mm



3812  
FLOHR  
3mm



3813  
FLOHR  
2mm



For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



**3816**  
CRYER  
27L



**3817**  
CRYER  
28R



**3814**  
CRYER  
39L



**3815**  
CRYER  
40R



**3818 OUT**  
BEIN-E77R  
3mm



3mm



**3820 IN**  
BEIN-E77  
4mm



**3822**  
HYLIN 2S



**3823**  
HYLIN 1STR



**3824**  
HYLIN 3S



**3749L**  
BURNARD 1L



**3749STR**  
BURNARD 2STR



**3749IN**  
BURNARD 3IN



**3749OUT**  
BURNARD 4OUT



**3749R**  
BURNARD 5R

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



3826  
COLE MAN  
1STR



3825  
COLE MAN  
1P



3827  
COLE MAN  
2LS



3828  
COLE MAN  
2RS



3835  
BARRY 11L



3836  
BARRY  
12R



3837  
BARRY  
17L



3838  
BARRY  
17R



3839  
WINTER 14L



3840  
WINTER 14R



For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



ELEVATORS

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



3807  
WWJ-L



3808  
WWJ-STR



3809  
WWJ-R



3887  
JK-1 2.5mm



3888  
JK-2 3.5mm



3889  
JK-3 3.5mm



3890  
JK-4 5mm



3891  
JK-5 5mm



3892  
JK-6 IN 2.7mm



3893  
JK-7 OUT 2.7mm



3878  
HYLIN  
12cm



3878  
HYLIN  
14cm



3878  
HYLIN  
16cm

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



**3854**  
FLOHR L

**3855**  
FLOHR R



**3856**  
L 2.5mm



**3857**  
R 2.5mm



**3858**  
L 3.5mm



**3859**  
R 3.5mm



**3860**  
L 4.5mm



**3861**  
R 4.5mm



**3862**  
STR 2.5mm



**3863**  
STR 3.5mm



**3864**  
L 2.5mm



**3865**  
R 2.5mm



**3866**  
L 3mm



**3867**  
R 3mm

## ELEVATORS

For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



3868

3869

3870  
SMALL L3871  
SMALL R3872  
LARGE L3873  
LARGE R3874  
2.5mm3875  
3mm3876  
4mm3877  
5mm



For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.



3901  
HEIDBRINK  
HB-1



3902  
HEIDBRINK  
HB-2



3903  
HEIDBRINK  
HB-3



3904  
HEIDBRINK  
HB-4



3905  
HEIDBRINK  
HB-5



3907  
APICAL  
9 STR



3908  
APICAL  
9R



3909  
APICAL  
9L



3916  
HOWARD-1



3917  
HOWARD-2



3918  
HOWARD-3



3919  
HOWARD-4



3920  
HOWARD-5




3906  
HEIDBRINK  
13/14



3921  
HOWARD  
2/3

For removing residual root tips of tooth.

DAVIS TEASERS



For pulling, lifting up the periosteum and making the tooth loosen from the periodontal ligament for easy extraction.

BLUE TITANIUM WITH SERRATED



3769 ACCU LUX STRAIGHT



3770 ACCU LUX CURVED



3771 ACCU LUX SPADE



3772 ACCU LUX DISTAL



3773 ACCU LUX MESIAL



3774 ACCU LUX-LARGE DISTAL



3775 ACCU LUX-LARGE MESIAL

# Forceps I

## SURGICAL SET UP 13 INS (Recommended)

1. 1812 6EZ SURGERY KNIVES KOR1/2
2. 2001 6EZ SURGICAL CHISELS 13K/TG
3. 3509 6EZ PERIOSTEALS PR3
4. 3513 6EZ PERIOSTEALS ALLEN9
5. 3514 6EZ PERIOSTEALS 24G
6. 4107 6EZ SURGICAL CURETTES LUCAS87
7. 4411 RONGEURS MICRO FRIEDMAN 15' (14.5CM)
8. 4618 STAINLESS STEEL RULER
9. 5505TN T/C MICRO CRILE-WOOD STR 15CM
10. 5905 KELLEY SCI, 16CM
11. 6420 RETRACTORS MINNESOTA
12. 6501 PLIERS ADSON 12CM
13. 6905 SCALPEL HANDLE 5



Extraction Forceps  
Root Splinters



## EXTRACTION FORCEPS AMERICAN PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.

All Handles are Unique Design  
Shown actual size.



**4801AF**  
1AF  
UPPER ANTERIORS



**4801**  
1  
UPPER CENTRALS  
LATERALS & CUSPIDS



**4802**  
5  
LOWER THIRD MOLARS



**4803**  
6  
SPLITTING FORCEPS



**4804**  
10H  
UPPER MOLARS  
HOOK HANDLE



**4805**  
15  
LOWER MOLARS



**4806**  
16S  
LOWER MOLARS



**4807**  
16  
LOWER MOLARS  
COW-HORN BEAK



**4808**  
17  
LOWER MOLARS



**4809**  
18L  
LOWER MOLARS  
LEFT SIDE



# EXTRACTION FORCEPS AMERICAN PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.

All Handles are Unique Design  
Shown actual size.



**4810**  
18R  
UPPER MOLARS  
RIGHT SIDE



**4811**  
23  
LOWER MOLARS  
COW-HORN BEAK



**4812**  
32  
UPPER MOLARS & PREMOLARS  
STRONG BEAK



**4813**  
53L  
UPPER MOLARS  
LEFT SIDE



**4814**  
53R  
UPPER MOLARS  
RIGHT SIDE



**4815**  
UNIVERSAL  
FOR UPPER TEETH & ROOTS



**4816**  
65  
UPPER INCISORS  
& ROOTS



**4817**  
69  
FRAGMENTS & SMALL ROOTS  
UPPER & LOWER



**4818**  
10S  
UPPER MOLARS  
STRAIGHT-HANDLE



**4819**  
88R  
UPPER MOLARS  
RIGHT SIDE



**4820**  
88L  
UPPER MOLARS  
LEFT SIDE



**4821**  
99C  
UPPER PREMOLARS  
& INCISORS



## EXTRACTION FORCEPS AMERICAN PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.



**4822**  
101  
UNIVERSAL  
FOR PREMOLARS



**4823**  
150  
UPPER PREMOLARS  
INCISORS & ROOTS  
CANINES



**4823AF**  
150AF  
UPPER PREMOLARS  
INCISORS & ROOTS  
CANINES



**4823AS**  
150AS  
UPPER PREMOLARS  
INCISORS & ROOTS  
CANINES



**4823XAS**  
150XAS  
UPPER PREMOLARS  
INCISORS & ROOTS  
CANINES



**4824**  
150A  
UPPER PREMOLARS  
INCISORS & ROOTS  
CANINES



**4825**  
150S  
UPPER PREMOLARS  
INCISORS & ROOTS  
FOR CHILDREN



**4826**  
151  
LOWER PREMOLARS  
INCISORS & ROOTS



**4826AF**  
151AF  
LOWER PREMOLARS  
INCISORS & ROOTS  
CANINES



**4826AS**  
151AS  
LOWER PREMOLARS  
INCISORS & ROOTS  
CANINES



**4826XAS**  
151XAS  
LOWER PREMOLARS  
INCISORS & ROOTS  
CANINES



**4827**  
151A  
LOWER PREMOLARS  
INCISORS & ROOTS

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.



**4828**  
151S  
LOWER PREMOLARS  
INCISORS & ROOTS  
FOR CHILDREN



**4829**  
210H  
UPPER WISDOMS  
WITH HOOK HANDLE



**4830**  
210S  
UPPER WISDOMS  
WITH STRAIGHT HANDLE



**4831**  
222  
LOWER WISDOMS



**4831AF**  
222AF  
LOWER WISDOMS



**4832**  
286  
UPPER ROOTS  
PREMOLARS & INCISORS



EXTRACTION FORCEPS  
ENGLISH PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.

All Handles are Unique Design  
Shown actual size.



4901

1  
UPPER CENTRALS  
& CANINES



4902

2  
UPPER LATERALS  
& BICUSPIDS



4903

3  
UPPER INCISORS  
& CANINES



4904

4  
LOWER INCISORS  
& CANINES



4905

7  
UPPER PREMOLARS



4906

8  
LOWER PREMOLARS



4907

13  
LOWER PREMOLARS



4908

17  
UPPER MOLARS  
RIGHT



4909

18  
UPPER MOLARS  
LEFT



4910

18A  
UPPER MOLARS



4911

19  
UPPER WISDOMS



4912

20  
LOWER WISDOMS

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.



**4913**  
21  
LOWER MOLARS



**4914**  
22  
LOWER MOLARS  
HAWKS BILL



**4919**  
30  
UPPER ROOTS  
& INCISORS



**4920**  
31  
LOWER ROOTS



**4921**  
32  
LOWER MOLARS



**4922**  
33  
LOWER ROOTS



**4925**  
33M  
LOWER ROOTS  
WITH LONG BEAK



**4927**  
44  
UPPER ROOTS



**4928**  
45  
LOWER ROOTS



**4929**  
51  
UPPER ROOTS



**4930**  
51A  
UPPER ROOTS  
NARROW BEAK



**4931**  
51C  
UPPER ROOTS



## EXTRACTION FORCEPS ENGLISH PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.



**4936**  
67  
UPPER WISDOMS



**4937AF**  
67AF  
UPPER WISDOMS



**4938**  
68  
LOWER ROOTS



**4939**  
73  
LOWER MOLARS  
(HAWKS BILL)



**4940**  
74  
LOWER INCISORS  
& ROOTS



**4941**  
74M  
LOWER INCISORS  
& ROOTS



**4942**  
74N  
LOWER INCISORS  
& ROOTS CANINES



**4942NF**  
74NF  
LOWER INCISORS  
& ROOTS



**4946**  
79  
LOWER WISDOMS



**4947**  
79A  
LOWER WISDOMS



**4950**  
86C  
LOWER MOLARS



**4951**  
87  
COWHORN BEAK

## EXTRACTION FORCEPS ENGLISH PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.



**4952**  
89  
UPPER MOLARS  
RIGHT



**4953**  
90  
UPPER MOLARS  
LEFT



**4954**  
94  
UPPER MOLARS  
RIGHT



**4955**  
95  
UPPER MOLARS  
LEFT



**4958**  
121  
UPPER WISDOMS



**4961**  
167  
UNIVERSAL PATER  
FOR LOWER TEETH  
MOLARS & PREMOLARS



**4962**  
168  
UPPER CANINES



**4963**  
159  
LOWER CANINES



**4964**  
233  
LOWER ROOTS  
PARALLEL BEAKS



**4965**  
235  
UPPER ROOTS  
PARALLEL BEAKS



**4966**  
MD-1  
UNIVERSAL  
UPPER WISDOM  
PREMOLARS & ROOTS



**4969**  
MD-4  
LOWER MOLARS



## EXTRACTION FORCEPS ENGLISH PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.



**4972**  
22 L  
LOWER WISDOMS  
LEFT

**4973**  
22 R  
LOWER WISDOMS  
RIGHT



**4974**  
F-1  
UPPER UNIVERSAL



**4975**  
F-2  
LOWER UNIVERSAL



**4976**  
F-3  
UPPER ANTERIOR



**4977**  
F-4  
LOWER ANTERIOR



**4977N**  
F-4N  
LOWER ANTERIOR



**4978**  
F-5  
LOWER ANTERIOR



**4979**  
F-6  
UPPER 3RD MOLAR



**4980**  
F-7  
LOWER 3RD MOLAR



**4981**  
F-8  
UPPER UNIVERSAL



**4982**  
F-9  
LOWER UNIVERSAL

EXTRACTION FORCEPS  
ENGLISH PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.



**4983**  
F-10  
LOWER MOLAR



**4984**  
F-11  
UPPER MOLAR  
LEFT



**4985**  
F-12  
UPPER MOLAR  
RIGHT



**4986**  
F-13  
UPPER MOLAR  
LEFT



**4987**  
F-14  
UPPER MOLAR  
RIGHT



**4988**  
F-15  
UPPER PREMOLARS  
UNIVERSAL



**4989**  
F-16  
UPPER ROOT



**4990**  
F-17  
LOWER ROOT



**4991**  
F-18  
UPPER ROOT  
FRAGMENT NARROW



**4991TN**  
F-18TN  
UPPER ROOT  
FRAGMENT NARROW



**4992**  
F-19  
LOWER ROOT  
FRAGMENT NARROW



**4992NW**  
F-19NW  
LOWER ROOT  
FRAGMENT NARROW

## EXTRACTION FORCEPS ENGLISH PATTERN

Forcep : For extracting teeth from alveolar bone.

American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter and Extraction Tweezer : For removing the residual part of tooth root.

Sponge and Dressing/Hemostats : For grasping off blood vessels and removing small pieces of root tips.



**4993**  
46L  
LOWER ROOT  
LONG BEAK



**4994**  
46XL  
LOWER ROOT  
EXTRA LONG



**4995**  
F-20  
UPPER PREMLAR  
1ST & 2RD MOLARS



**4996**  
49  
UPPER ROOT



**4997**  
51XL  
UPPER ROOT



ALL HANDLES ARE UNIQUE DESIGN



CHILDREN PATTERNS



**5001**  
158 UPPER MOLARS  
UNIVERSAL



**5002**  
159 UPPER TEETH  
UNIVERSAL



**5003**  
160 LOWER MOLARS  
UNIVERSAL



**5004**  
161 LOWER CENTRALS &  
ROOTS



**5005**  
162 UPPER CENTRALS



**5006**  
163 UPPER CENTRALS



**5007**  
37 UPPER INCISORS &  
CANINES



**5008**  
39 UPPER MOLARS  
UNIVERSAL



**5009**  
123 UPPER TEETH  
UNIVERSAL



**5010**  
138 UPPER UNIVERSAL &  
ROOTS



**5011**  
157 UPPER MOLARS  
UNIVERSAL

PEDODONTIC PATTERNS



**5012**  
3 UPPER MOLARS



**5013**  
5 LOWER INCISORS



**5014**  
6 LOWER MOLARS



**5015**  
7 LOWER ROOTS



**5016**  
137 UPPER INCISORS &  
CANINES



**5017**  
139 UPPER PREMOLARS



**5018**  
51S UPPER ROOTS



**5019**  
150-SK UPPER UNIVERSAL



**5020**  
151-SK LOWER UNIVERSAL



**5022**  
B23 LOWER MOLARS



For extracting root tips of the tooth where normal forceps cannot reach.



5101  
301  
LOWER ROOTS



5102  
300  
UPPER ROOTS



5103  
ALBRECHT 16cm



5103  
ALBRECHT 16cm



5104  
ALBRECHT 16cm



# Crown Instruments J

## CROWN AND BRIDGE SET UP 19 INS (Recommended)

1. 1006 3EZ EXPLORERS EXP 5
2. 1104 1EZ COLOR PROBE CP12
3. 1310 6EZ SCALERS/CURETTES C11, CEMENT REMOVER
4. 1318 6EZ SCALERS/CURETTES MC1, CEMENT REMOVER
5. 2102 4EZ HATCHETS H13/14
6. 2103 4EZ HATCHETS H15/16
7. 2206 4EZ EXCAVATORS 17
8. 2501 4EZ CARVERS IPC
9. 2702 4EZ CARVING KNIVES 21
10. 2807P 4EZ ACCURATOR-PLACER
11. 2910 4EZ PLASTIC FILLING TUFTS4
12. 3109 4EZ CORD PACKER CSI1
13. 4602T CROWN SCISSOR CVD T/C
14. 4706 FORCEP PAPER
15. 6517 FORCEP LONDON-COLLEGE
16. 6807 SPATULAS
17. 6911 6EZ MIRROR HANDLE #6 C/S
18. 6951 #5 MIRROR C/S FRONT SURFACE
19. 7502N CASSETTE 16 HOLD DOUBLE DECK
20. N CEMENT REMOVER, 1310(C11) + 1318(1)



Crown Instruments  
Crown Removers  
Trays

Crown Scissors – For removing crowns permanently by cutting off the seal between tooth and crown after splitting with a bur.

Crown Spreader – For splitting and spreading the crown open for removal.

Crown Remover – For gripping onto the crown and removing from the tooth



4601  
STRAIGHT  
SAW EDGE  
10.5cm  
4601T  
T/C



4602  
CURVED  
SAW EDGE  
10.5cm  
4602T  
T/C



4603  
UNIVERSAL  
12cm



4604  
IWANSON



4605  
CROWN SPREADER



4606  
BOLEY 14cm



4607  
CROWN REMOVER  
PLIERS



4608  
CROWN SPREADER  
FORCEPS



4610  
CROWN GRIPPER



For removing crowns permanently by cutting off the seal between tooth and crown after splitting with a bur.  
 Measuring Instruments – For taking accurate measurements in specific places around the bone.



**4613**  
 CASTROVIEJO  
 STR 8.5cm



**4613L**  
 CASTROVIEJO  
 STR 17cm



**4613LA**  
 CASTROVIEJO  
 CVD 17cm



**4613RA**  
 CASROVIEJO  
 17cm



**4614**  
 RIDGE MAPPING BONE width 0 ~25mm  
 14cm



**4615**  
 CROWN STR



**4615T**  
 CROWN STR T/C



**4616**  
 BEERENDONT



**4616D**  
 DOUBLE CALIPER



**4617**  
 CROWN GRIPPER



**4617N**  
 CROWN GRIPPER



**4618**  
 RULER  
 15cm / 6in



**4625**  
 CALIPER DOUBLE



**4627**  
 ENDO FINGER ROLLER





4701  
CROWN REMOVERS



4702  
CROWN REMOVERS



4703  
CROWN REMOVERS



4704  
CROWN REMOVERS



4705  
CROWN SPREADER



4706  
ARTICULATING PAPER  
FORCEPS 15cm



4707  
MILLER 15cm



4708 PLASTER CUTTING PLIER 17cm

4709 PLASTER CUTTING PLIER 20cm



4710  
CROWN SPREADER  
FLAT



4710T  
CROWN SPREADER  
TEETH

For forming a crown frame.



SOLID PARTIAL  
SET OF 10

7301  
UPPER

7301  
LOWER



PERFORATED PARTIAL  
SET OF 10

7302  
UPPER

7302  
LOWER

	UPPER	U-1(XS)	U-2(S)	U-3(M)	U-4(L)	U-5(XL)
	A	65	68	71	73	78
	B	52	54	58	59	63

	LOWER	L-1(XS)	L-2(S)	L-3(M)	L-4(L)	L-5(XL)
	A	66	70	73	76	78
	B	50	52	56	58	60

----- 7301 -----  
----- 7302 -----



SOLID FULL DENTURE  
SET OF 8

7303  
UPPER

7303  
LOWER



PERFORATED FULL DENTURE  
SET OF 8

7304  
UPPER

7304  
LOWER

	UPPER	U-1(S)	U-2(M)	U-3(L)	U-4(XL)
	A	65	66	72	77
	B	58	63	61	68

	LOWER	L-1(S)	L-2(M)	L-3(L)	L-4(XL)
	A	65	72	75	83
	B	56	60	61	62

----- 7303 -----  
----- 7304 -----

TRAYS FOR CROWN & BRIDGE



7307  
UPPER

7307  
LOWER



7308  
UPPER

7308  
LOWER



7309  
UPPER RIGHT  
LOWER LEFT

7309  
UPPER LEFT  
LOWER RIGHT



7310  
ADJUSTABLE  
PERFORATED



7311  
ADJUSTABLE  
SOLID



7312  
BITE CHECK TRAY

\* All Dimension are in 'mm'

288 X 187 X 29mm(WITH LID)



7401  
SOLID BASE  
WITHOUT SEPARATOR

200 X 110 X 15mm



7407  
SCALERS TRAY



7408  
6" X 10"



7410  
215 X 127 X 13mm  
215 x 165 x 16mm



7412  
ENDODONTIC BOX



7423  
SHARPENING STONE  
INDIA SLIP



7424  
SHARPENING STONE  
ARKANSAS FLAT



7425  
TRAY IMPLANT



7426  
TRAY IMPLANT



# Orthodontic **K**

## ORTHODONTIC SET UP 12 INS (Recommended)

1. 1203 6EZ SCALERS/CURETTES H5/33
2. 5603 WIRE CUTTER 15 T/C
3. 5607 WIRE CUTTER STR THIN SOFT NON T/C
4. 5611 PLIER WEINGART THINNER
5. 5617 PLIER ADHESIVE REMOVE
6. 5718 BIRD BEAK 12.5CM WITH SOFT WIRE CUTTER
7. 5620A BRACKET REMOVER ANGLE
8. 5723 ORTHODONTIC PLIER ADERER 12CM
9. 5805 ORTHODONTIC PUSHER 14CM
10. 6911 6EZ MIRROR HANDLE #6 C/S
11. 6951 #5 MIRROR C/S FRONT SURFACE
12. 7502N CASSETTES 16 HOLD DOUBLE DECK



Cutters & Pliers  
Matrix Retainers



For cutting wire during orthodontic procedures.



5601, 5601S, 5601L  
Universal Cutter  
& Distal End With  
T/C



5602



5602H  
Distal End Cutter  
& Hold



5603, 5608, 5609S  
Cutter 15°



5625  
15°

For cutting or bending wire during orthodontic procedures.

**BAND REMOVER:** A special plier that is used to remove orthodontic bands from teeth.

**BRACKET REMOVER:** For use in the removal of the bracket from the tooth.

**BAND SEATER:** For use in placement of the band around the tooth.



5610  
WEINGART  
(THICKER)



5611  
WEINGART  
(THICKER)



5612  
HOW  
(STRAIGHT)



5613  
HOW  
(STRAIGHT & LONG)



5614  
HOW  
(ANGLED)



5615  
HOW  
(ANGLED)



5616  
PLIER  
(BAND SEATING)



5617  
PLIER  
(ADHESIVE REMOVER)



5618  
PLIER  
(BAND REMOVER)



5620  
PLIER  
(BRACKET REMOVER)



5620A  
PLIER  
BRACKET REMOVER



5621  
PLIER



5622  
HOW 15cm



5623  
NERVE CANAL NARROW  
12.5 cm



5624  
HARD WIRE CUTTER  
STRAIGHT



5625  
HARD WIRE CUTTER  
CVD 15°



For bending wire during orthodontic procedures.



**5701**  
WALDSACHE  
15cm



**5702**  
UNIVERSAL  
16cm



**5703**  
NANCE  
13.5cm



**5704**  
ADERER  
14.5



**5705**  
GOSLEE  
ONE SERRATED 14cm



**5706**  
GOSLEE  
DOUBLE SERRATED 14cm



**5709**  
ARROW CLASP  
BENDING 12cm



**5710**  
ARROW CLASP  
FORMING 12cm



**5712**  
PEESO  
13.5cm



**5714**  
POINTED JAWS  
14cm



**5716**  
SLIM NOSE JARABEK  
14cm



**5717**  
13cm



**5718**  
BIRD BEAK  
12.5cm



**5718L**  
BIRD BEAK  
15.5cm



**5721**  
SCHWARZ  
13.5cm



**5723**  
THREE JAW  
ADERER 12cm



**5723N**  
THREE JAW  
ADERER SMALL 11.5cm



**5724**  
TWEED  
12cm



**5725**  
O BRIEN  
13.5cm



**5726**  
WIRE BENDING  
12.5cm



**5727**  
YOUNG  
12.5cm



**5728**  
CLASP FORMING  
14cm



**5729**  
ADAMS  
12.5cm



**5730**  
FREVERT  
12.5cm

For bending wire during orthodontic procedures.



5731  
NANCE  
12.5cm



5732  
GROVES BENDING  
12.5cm



5733  
BENDING PLIER 0.75mm  
12.5cm



5734  
BENDING PLIER 1mm  
12.5cm



5741  
ORTH RACK 22cm





Commonly used instruments in orthodontic procedures.



5801N  
ORTH



5801  
BAND SEATING &  
REMOVING



5802



5803



5804  
ELASTIC



5805  
PUSHER  
MERSHON BAND  
14cm



5806  
BOND SETTER  
15cm



5807  
13.5cm



5808  
14cm



5809  
SEPARATING PLIERS  
15.5cm



5810  
MATHIEU PLIERS  
14cm



5812  
MATHIEU LIGATURE  
PLACING PLIERS  
14cm



5814  
MOSQUITO PLIERS  
12.5cm



5816  
BRACKET PLACING  
12.5cm



5817  
BRACKET  
POSITIONING GAGE



5818  
LIGATURE TYING  
PLIERS



5819  
BRACKET HOLDING  
15cm



5820  
BOND SEATING  
13cm



5821  
MOSQUITO FORCEPS  
12.5cm

Commonly used instruments in orthodontic procedures.



5830  
LDP  
LIGATURE DIRECTOR/PLUGGER



5831  
LOT  
LIGATURE DIRECTOR/TUCKER



5832  
TWEezer



5833  
SCALER/LIGATURE DIRECTOR



5840  
U15/F  
SCALER/BAND PUSHER



5841  
PICKER/LIGATURE DIRECTOR



5842  
TWISTER

MATRIX RETAINERS

For filling in the amalgam materials when cavity occurs at the side of tooth.

Matrix Retainer – To hold metal bands or strips around the tooth creating a temporary mold while filling material is packed into place.



6201  
IVORY 8  
5mm

6203  
IVORY 8  
7mm



6205  
IVORY 8

6204  
IVORY 8N



6206  
SIQVELAND  
NORROW

6207  
SIQVELAND  
WIDE



6208  
1  
TOFFLEMIRE  
UNIVERSAL



6209  
3  
TOFFLEMIRE  
SENIOR



6210  
2  
TOFFLEMIRE  
JUNIOR



6212  
IVORY



6213  
IVORY



6214  
IVORY  
SEPARATOR



6216  
PAPER HOLDING FORCEP  
DOUBLE JAW







The ergonomic hold - The design of the tissue and cotton pliers prevents slipping while reducing the weight by 20%.

UNIQUE GRIP DESIGN



UNIQUE GRIP DESIGN



6504  
SEMKEN-TAYLOR 12.5cm



6505  
SEMKEN-TAYLOR 12.5cm



6506  
SEMKEN-TAYLOR 12.5cm



6507  
12cm



6508  
12cm



6509  
12cm



6510  
12cm

The ergonomic hold - The design of the tissue and cotton pliers prevents slipping while reducing the weight by 20%.



6511  
SUTURE 16cm



6511N  
SUTURE 16cm



6512  
GOLDMAN-FOX  
(CRANE-KAPLANE)  
15cm



6513  
GOLDMAN-FOX  
(CRANE-KAPLANE)  
15cm



6515  
SELF-LOCKING 15cm



6517  
LONDON-COLLEGE 15cm



6517L  
LONDON-COLLEGE  
(SELF-LOCKING)  
15cm



6518  
MERIAM 15cm



6519  
PERRY 13cm



6520  
LONDON-COLLEGE  
(SELF-LOCKING)  
15cm



6521  
YASARGIL 18cm



6522  
SELF-LOCKING 16cm



6523  
SELF-LOCKING



6525  
TISSUE FORCEP T/C  
GERALD 18cm



6526  
TOMOR 16cm



6526L  
TOMOR 20cm





The ergonomic hold - The design of the tissue and cotton pliers prevents slipping while reducing the weight by 20%.



6527  
TF16 15cm



6528  
TF50L  
SELF-LOCKING 15cm



6531 DEBAKEY ATRAUMATIC 18cm  
6531C DEBAKEY ATRAUMATIC CVD 18cm



6529  
ALLISON BABY  
12.5cm



6530  
ALLISON 15cm

For laboratory use for carving wax.



For laboratory use for carving wax.



TOWA INSTRUMENTS



# SUCTION TUBES

Suction system tubes for removing saliva, blood, pus, and solids that accumulate during dental procedures.



6714 2mm

6715 2.5mm

6716 3mm

6717 4mm

6735 5mm



6731  
SUCTION TUBE 1.5mm



6732  
SUCTION TUBE 2.5mm



6733  
SUCTION TUBE 3mm



6734  
SUCTION TUBE 4mm





Towel Clamps – Forceps used to secure drapes to the skin and to hold tissue.



7101  
BACKHAUS  
8cm



7102 10.5cm  
7103 9cm  
7104 11cm



7105  
BACKHAUS  
13.5cm



7107  
16cm



7108  
CHEATTLE 16cm



7109  
20cm



7110  
20cm/8"  
30cm/12"



7111  
ROGGE  
23.5cm

Instruments Not Included.



**7501**  
**5 INSTRUMENT CASSETTE**  
 (2.5" X 7" X 0.75")  
 (64mm X 178mm X 18mm)



**7502 10 INS**    **7502N 16 INS**  
**8/10/16 INSTRUMENT CASSETTE**  
 (4.75" X 8" X 1.25")  
 (120mm X 203mm X 32mm)



**7503 6 INS**    **7503N 12 INS**  
**8/12 INSTRUMENT CASSETTE**  
 (3.75" X 8" X 1.25")  
 (95mm X 203mm X 32mm)



**7505**  
**20 INSTRUMENT CASSETTE**  
 (8" X 11" X 1.25")  
 (203mm X 280mm X 32mm)



**7505-1**  
**SCISSOR CASSETTE CLIP**



**7505-2**  
**SUCTION CASSETTE CLIP**



Internal and external flow scaler inserts.



7601 25KHZ UNIVERSAL No 10/50 EXTERNAL FLOW METAL GRIP



7602 25KHZ SLIM LINE No 100 EXTERNAL FLOW METAL GRIP



7603 30KHZ UNIVERSAL No 10/50 EXTERNAL FLOW METAL GRIP



7604 30KHZ SLIM LINE No 100 EXTERNAL FLOW METAL GRIP



7606 25KHZ STREAM LINE No 10 INTERNAL FLOW PLASTIC GRIP



7607 25KHZ SLIM LINE No 100 INTERNAL FLOW PLASTIC GRIP



7610 30KHZ STREAM LINE No 10 INTERNAL FLOW PLASTIC GRIP



7611 30KHZ SLIM LINE No 100 INTERNAL FLOW PLASTIC GRIP



7605 25KHZ BEAVERTAIL No 3 INTERNAL FLOW PLASTIC GRIP



7609 30KHZ BEAVERTAIL No 3 INTERNAL FLOW PLASTIC GRIP